

## MICROPLASTIC DETECTION

### ECSIN LAB RESEARCH AND METHODS VALIDATION

The vast array of different plastic types available on today's consumer market makes the qualitative or quantitative analysis of microplastics extremely challenging and there are no officially recognized methods available by now.

However thanks to **ECSIN LAB progress in scientific research, specific methods applicable to different matrix types** - including detergents, cosmetics, drinking water, milk, beverages, mineral salts, fish products and environmental matrices (e.g. waste water, soil, sludge) - were **developed and validated** through:

- Verification of process and possible contamination through blank analysis
- Recovery assessment of the number of microplastics added
- Assessment of the integrity (size and shape) of added microplastics
- Verification of chemical nature of added microplastics

**Federico Benetti, PhD**  
ECSIN Division Director  
f.benetti@ecamricert.com

**MICROPLASTIC DETECTION**  
ECSIN LAB **ACCREDITATION**



## VALIDATED PROCEDURES FOR

Salt  
Drinking Water  
Milk  
Soft drinks  
Mollusks  
Fish



Waste water  
Fresh water / Sea water  
Sludge  
Soils  
Digestate



Detergents  
Cosmetics



**ecsinn**  
european center  
for the sustainable impact  
of nanotechnology

# MICROPLASTIC DETECTION

## ECSIN LAB RESEARCH AND METHODS VALIDATION

### Sample preparation STEP 1

The first step is the degradation of the organic matrix in order to isolate microplastics

LIQUID SAMPLES  
SOLID SAMPLES

CHEMICAL AND/OR  
ENZYMATIC DIGESTION



### ARTIFACT PREVENTION

- Alteration or particle loss
- External contamination

### Filtration STEP 2

Microplastics have to be completely recovered on filter and detectable.

FILTRATION



- Filtration in cleanroom ISO 7 class
- Plastic-free items
- Appropriate cleaning procedure before filtering

### μFTIR / chemical imaging STEP 3


#### μFTIR SAMPLING MODES

- Reflectance
- Transmittance
- ATR (attenuated total reflection)

#### μFTIR / CHEMICAL IMAGING

Analytical capability to create a **visual image of components distribution** from simultaneous chemical identification and particle size detection:


- Generate **chemical images** showing differences
- Each **pixel corresponds** to IR spectrum



# MICROPLASTIC DETECTION

## FILTRATION STEP

Matrix degradation without affecting **MICROPLASTICS**, both in size and composition.



SIZE DISTRIBUTION PARAMETERS	
Number of particles	206
Minimum diameter (µm)	89
First quartile (µm)	97
Median (µm)	99
MAD (µm)	4
Average (µm)	99
Standard deviation (µm)	4
Third quartile (µm)	101
Maximum diameter (µm)	115

**MICROPLASTICS** have to be completely recovered on filter.

		SAMPLE 1	SAMPLE 2	SAMPLE 3	SAMPLE 4	SAMPLE 5	SAMPLE 6
Polyethylene (10-106 µM)	Number of spiked mps in sample	335	834	355	416	671	631
	Recovery in number	289	844	408	446	700	587
	RECOVERY %	86%	101%	115%	107%	104%	93%




## Case study

# MICROPLASTICS DETECTION IN CRUSTACEANS




**Crustaceans: flagship species for deep-sea environments.**

Analyzed **crustaceans** present microplastics particles in the size range **36 – 2556  $\mu\text{m}$** . Particles are mainly composed of **polyethylene (PE)** and **polipropilene (PP)**. A second analysis examined particles smaller than 150  $\mu\text{m}$  that identify microplastics in the size range of in the size range 36 – 132  $\mu\text{m}$ .



**CRUSTACEANS: 12 particles / g**

SIZE RANGE: 36 – 2556  $\mu\text{m}$



**CRUSTACEANS: 10 particles / g**

SIZE RANGE: 36 – 132  $\mu\text{m}$




## Case study

### MICROPLASTICS DETECTION IN TOMATO SUPPLY CHAIN CONTAMINATION



Analyzed tomato present microplastics particles in the size range 33 – 200  $\mu\text{m}$ . Particles are mainly composed of **polyethylene terephthalate (PET)**. A second analysis examined particles smaller than 150  $\mu\text{m}$  that identify microplastics in the size range of in the size range 33 – 130  $\mu\text{m}$ .



## Case study

### MICROPLASTICS DETECTION IN SOAP CUBE POST-PROCESS CONTAMINATION



Analyzed soap cube present microplastics particles in the size range 37 – 139  $\mu\text{m}$ . Particles are mainly composed of **polyethylene (PE)**.

