

Bacterial 16S metagenomics by high-throughput sequencing: a tool for the etiological diagnosis of pneumonia in cattle?

- Study done by Fabien Gregoire, ARSIA
- Reported by Christian Quinet FESASS

In 2020 :

1490 bovine necropsied (all ages)

- ❖ 1177 bovine necropsied (0-6 months)
- ❖ In this age group, 293 animals with respiratory lesions

Approximately 80% of respiratory cases elucidated

Based on:

- ✓ Macroscopic lesions
- ✓ Histology
- ✓ Culture (Gassner - blood agar - chocolate - CNA - PAM)
- ✓ PCR (6 pathogens)




FESASS

Respiratory diseases at ARSIA: some numbers




PCR results


results on autopsy lungs (2018-2019-2020):



Results of autopsy lung cultures (2020):





Aerobic culture



Mycoplasma culture

16 S Sequencing Technology



Kazumasa Fukuda, 2016. DOI:[10.7888/juoeh.38.223](https://doi.org/10.7888/juoeh.38.223)

Lung sample

Storage -20°C

Ethanol disinfection (min 24h)

PBS rinsing and heart sampling

DNA extraction


DNA quantification (Quantus™ Promega)


PCR 16S regions V4 and V7 (ThermoFisher primers)

Illumina Miseq sequencing

Analysis of results with Geneious software
(reference base: NCBI)

Processing of results





**Calves 1-6 month
From necropsy
room**



**With pneumonia lesions
N = 174**


**58 « resolved »
116 « not resolved »**

**Healthy lungs
N = 24**




<i>Trueperella pyogenes</i> strain ATCC 19411	16S ribosomal RNA	30321	18.73%
<i>Fusobacterium necrophorum</i> strain JCM 3718	16S ribosomal RNA	29902	18.47%
<i>Bacteroides heparinolyticus</i> strain HEP	16S ribosomal RNA	22274	13.76%
<i>Porphyromonas endodontalis</i> strain JCM 8526	16S ribosomal RNA	18678	11.53%
<i>Mycoplasma bovirhinis</i> strain NBRC 14857	16S ribosomal RNA	12178	7.52%
Unclassified (conflicting taxonomy)		8538	5.27%
<i>Porphyromonas somerae</i> DSM 23386 strain JCM 13867	16S ribosomal RNA	7945	4.91%
<i>Porphyromonas crevioricanis</i> strain NUM 402	16S ribosomal RNA	6889	4.25%
<i>Helcococcus ovis</i> strain s840-96-2	16S ribosomal RNA	5002	3.09%
<i>Parvimonas micra</i> strain JCM 12970	16S ribosomal RNA	4696	2.90%
<i>Peptoniphilus indolicus</i> ATCC 29427	16S ribosomal RNA	3407	2.10%
<i>Mycoplasma dispar</i> strain 462/2	16S ribosomal RNA	3005	1.86%
<i>Filifactor alocis</i> strain ATCC 35896	16S ribosomal RNA	2509	1.55%
Unclassified (no match)		1603	0.99%
<i>Bacteroides pyogenes</i> strain DSM 20611	16S ribosomal RNA	1245	0.77%
....			
....			
....			

Case/control study



Case/control study





FESASS

Case/control study



Pathogenic bacteria routinely tested :

Mannheimia haemolytica

Mycoplasma bovis

Pasteurella multocida

Trueperella pyogenes

Histophilus somni

Pathogenic bacteria not routinely tested :

*Mycoplasma dispar**

*Helcococcus ovis**

*Ureaplasma diversum**

*Mycoplasma bovirhinis**

Pathogenicity/identification not confirmed :

Streptobacillus hongkongensis

*Peptinophilus indolicus**

*Bacteriodes pyogenes**

Parvimonas micra

Porphyromonas somerae

Pseudomonas trivialis

Trueperella bernadiae (= *Trueperella pyogenes*?)

...

HELCOCOCCUS OVIS ISOLATED FROM SUPPURATIVE ARTHRITIS IN A LAMB

Lucienne Downs, District Veterinarian, Central Tablelands Local Land Services

Posted Flock & Herd August 2020

Helcococcus ovis, an Emerging Pathogen in Bovine Valvular Endocarditis

September 2008 · *Journal of Clinical Microbiology* 46(10):3291-5 · [Follow journal](#)

DOI: [10.1128/JCM.00867-08](https://doi.org/10.1128/JCM.00867-08)

Front. Vet. Sci., 11 February 2022 | <https://doi.org/10.3389/fvets.2021.756438>



Biological Characteristics and Pathogenicity of *Helcococcus ovis* Isolated From Clinical Bovine Mastitis in a Chinese Dairy Herd

Kai Liu¹, Zhaoju Deng², Limei Zhang¹, Xiaolong Gu¹, Gang Liu¹, Yang Liu¹, Peng Chen¹, Jian Gao³, Bo Han⁴ and Weijie Qu¹

Association of *Ureaplasma diversum* with reproductive disease in cattle

JM Diaz [✉](#), A Prieto, G López, P Díaz, C López, LÁ Quintela, ...show all

Pages 249-256 | Received 08 Oct 2018, Accepted 18 May 2019, Accepted author version posted online: 26 May 2019, Published online: 03 Jul 2019

[Download citation](#) <https://doi.org/10.1080/00480169.2019.1623733> [Check for updates](#)

* Taxons confirmed by PCR (ULiège) or bacteriology

Resolving these 'unresolved' cases by 16S?

CI-19-164888	CI-19-164888-001 AT Avis du pathologiste sur ce dossier Commentaire: CI-19-164888-001 AT Résultat de l'autopsie OK CI-19-164888-001 BA Culture Aérobic Flore non spécifique	Helcococcus ovis strain s840-96-2 16S ribosomal RNA 11090 28.53%	Unclassified (no match) 9166 23.58%	Trueperella bernardiae strain LCDC 89/0504 16S ribosomal RNA 7852 20.20%	Ignatzschineria indica strain FFA1 16S ribosomal RNA 3087 7.94%	Streptobacillus hongkongensis strain HKU33 16S ribosomal RNA 3017 7.76%	Pseudomonas silesiensis strain A3 16S ribosomal RNA 842 2.17%	Streptobacillus felis strain 131000547 16S ribosomal RNA 676 1.74%	Ureaplasma diversum strain A417 16S ribosomal RNA 548 1.41%	Unclassified (conflicting taxonomy) 461 1.19%	Moraxella cuniculi strain CCUG 2154 16S ribosomal RNA 454 1.17%	...
CI-19-167815	CI-19-167815-001 AT Avis du pathologiste sur ce dossier Commentaire: CI-19-167815-001 AT Résultat de l'autopsie OK CI-19-167815-001 BA Culture Aérobic Flore non spécifique	Ureaplasma diversum strain A417 16S ribosomal RNA 60782 71.88%	Pseudomonas chlororaphis subsp. aurantiaca strain NCIB 10068 16S ribosomal RNA 9178 10.85%	Aeromonas salmonicida subsp. masoucida strain NBRC 13784 16S ribosomal RNA 4642 5.49%	Unclassified (no match) 3560 4.21%	Unclassified (conflicting taxonomy) 2962 3.50%	Pseudomonas tolaasii strain ATCC 33618 16S ribosomal RNA 863 1.02%	Erwinia persicina strain NBRC 102418 16S ribosomal RNA 510 0.60%	Pseudomonas stutzeri ATCC 17588 = LMG 11199 16S ribosomal RNA 492 0.58%	Mycoplasma dispar strain 462/2 16S ribosomal RNA 386 0.46%	Serratia liquefaciens strain ATCC 27592 16S ribosomal RNA 161 0.19%	...
CI-19-231163	CI-19-216040-003 BO Biothèque oui CI-19-231163-001 AT Avis du pathologiste sur ce dossier Commentaire: CI-19-231163-001 AT Résultat de l'autopsie OK CI-19-231163-001 BA Culture Aérobic Envahissement par Proteus sp	Pseudomonas mendocina strain ATCC 25411 16S 15433 29.17%	Mycoplasma bovirhinis strain NBRC 14857 16S 9950 18.80%	Pseudomonas aeruginosa strain DSM 50071 16S 7416 14.01%	Unclassified (conflicting taxonomy) 7272 13.74%	Bacteroides pyogenes strain DSM 20611 16S 5992 11.32%	Xenorhabdus hominickii strain KE01 16S 1557 2.94%	Proteus mirabilis strain ATCC 29906 16S 1276 2.41%	Proteus alimentorum strain 08MAS0041 16S 1090 2.06%	Mycoplasma dispar strain 462/2 16S 846 1.60%	Trueperella pyogenes strain ATCC 19411 16S 593 1.12%	...
CI-20-056519	CI-20-056519-001 AT Avis du pathologiste sur ce dossier Commentaire: CI-20-056519-001 AT Résultat de l'autopsie OK CI-20-056519-001 BA Culture Aérobic Contaminants CI-20-056519-001 BA Culture Aérobic Flore non spécifique	Mannheimia haemolytica strain NCTC 9380 16S ribosomal RNA 18136 52.93%	Unclassified (no match) 4692 13.69%	Pseudomonas silesiensis strain A3 16S ribosomal RNA 3545 10.35%	Pseudomonas tolaasii strain ATCC 33618 16S ribosomal RNA 1322 3.86%	Serratia liquefaciens strain ATCC 27592 16S ribosomal RNA 996 2.91%	Gardnerella vaginalis ATCC 14018 = JCM 11026 16S ribosomal RNA 677 1.98%	Unclassified (identity too low) 516 1.51%	Escherichia marmotae strain HT073016 16S ribosomal RNA 490 1.43%	Unclassified (conflicting taxonomy) 384 1.12%	Corynebacterium coyleae strain DSM 44184 16S ribosomal RNA 230 0.67%	...
CI-21-012910	CI-21-012910-001 AT Avis du pathologiste sur ce dossier Commentaire: CI-21-012910-001 AT Résultat de l'autopsie OK	Unclassified (conflicting taxonomy) 16845 39.27%	Unclassified (no match) 6093 14.21%	Pseudomonas tolaasii strain ATCC 33618 16S ribosomal RNA 4910 11.45%	Staphylococcus capitis strain JCM 2420 16S ribosomal RNA 4093 9.54%	Ureaplasma diversum strain A417 16S ribosomal RNA 2981 6.95%	Cutibacterium acnes strain ATCC 6919 16S ribosomal RNA 1196 2.79%	Staphylococcus aureus strain ATCC 12600 16S ribosomal RNA 935 2.18%	Staphylococcus pasteurii strain ATCC 51129 16S ribosomal RNA 686 1.60%	Corynebacterium tuberculostrictum strain Medalle X 16S ribosomal RNA 502 1.17%	Polygonibacillus indicireducens strain In2-9 16S ribosomal RNA 457 1.07%	...

✓ Approximately 50% of "unsolved" cases solved with 16S



FESASS



Resolving these 'unresolved' cases by 16S?

Limitations of the 16S method:

- Identification of taxa
- Lot of different bacteria
- Considering as pathogenic or not pathogenic
- Quality of samples and "invasion" by contaminating bacteria
- What about the viral sphere?

Conclusions et perspectives



Some (nice) surprises



Some (nice) challenges



Some nice surprises!

- ✓ Confirmation of pathogens researched in routine
- ✓ Identification of "new" pathogens of interest
- ✓ Opening new areas for exploration
- ✓ Usable in routine?
 - Currently complex and time-consuming at the technical and data processing levels --> bioinformatics flow to be consolidated
 - Only targets the bacterial sphere
 - Interesting for "difficult" cases

Most interest is the exploration of underestimated or unknown pathogens => Identification of new/not yet known risks for respiratory diseases in bovine

Some nice challenge!

- F
New case/control study for Helcococcus ovis,
Mycoplasma dispar et Ureaplasma diversum :
 - Prevalence (with homemade PCR)
 - Epidemiological study
 - Standardised histology
- C
- E
of
interest



Some nice challenge!

- Further analysis of the results :
 - Quantitative aspect
 - Co-infections
 - Comparison of reference databases
- Confirmation of taxa of interest
- Epidemiologic and pathological study on confirmed taxas of interest
- Investigation of other systems/pathologies (abortions, metritis, ...)
- Investigation of the viral sphere



A team work !

ARSIA

Molecular Biology Department : Isabelle, Pascale et Julie

Pathology Department: Thierry, Frédéric and staff

Julien Evrard, coordinator of all Arsia projects

Fabien Gregoire, manager of this project

ULiège, Veterinary University of liège

Prof. MM Garigliani

Thanks for attention !